



4xM3x0.5 - 6H

Ø26,20

4xM3x0.5 - 6H

4xM3x0.5 - 6H

D-Sub 15pins male

-All dimension are in mm.  
 -All Tolerances are accordance with ISO 2768-1 class f, unless otherwise specified.  
**Tolerances:**  
 -Internal corners R0.2 max.  
 -Break edges 0,1 - 0,25  
 -√Rz 0,4 max  
 -Remove all sharp edges and burrs

Date.	Rev.	QTY.
15-06-2022	0	1
Designer.	Mass.	
JJ X-Ray A/S	0,993 kg	
Surface Treatment.		
Material.		

**JJ X-RAY**  
 Danish Science Design

Project Name  
 JJ X-Ray ATC8HV Slit

Description.

Comments.	
Scale:	Drawing Size.
⊕	A3
CAD Drawing Number	
22033000	