



M6x1 - 6H

M6x1 - 6H

M6x1 - 6H

M6x1 - 6H

-All dimension are in mm.
 -All Tolerances are accordance with ISO 2768-1 class f, unless otherwise specified.
Tolerances:
 -Internal corners R0.2 max.
 -Break edges 0,1 - 0,25
 -√Rz 0,4 max
 -Remove all sharp edges and burrs

Date: 06-02-2019
 Designer: JJ X-Ray A/S
 Surface Treatment:

Rev. 0
 Mass: 18,583 kg
 QTY.

JJ X-RAY
 Danish Science Design
 Project Name: IB C160 Air slit.
 Description: IB C160 Air slit

Comments:
 Scale: 1:1
 Drawing Size: A2

CAD Drawing Number: 12053002